

Get Your High Quality Instruments


tkmed
INSTRUMENTS



DENTAL INSTRUMENTS

TKMED Instruments is a brand for manufacturing quality instruments. An eminent name among Disposable & Reusable Dental instruments manufacturers in Pakistan. The **TKMED Instruments** features a wide range of products details are available on website and catalogue.

www.tkmed-instruments.com



01

Extracting Forceps


Page 01-31



02

Root Elevators


Page 32-50



03

Filling Instruments


Page 51-55



04

Dental Probes


Page 56-59



05

Periodontia Instruments
Mallets

Page 60-67



06

Curettes
Scalers

Page 68-74

TABLE OF CONTENTS

07



Page 75-77

Modelling Instruments
Measuring Instruments


08



Page 78-79

Wax / Plaster Knives


09



Page 80-85

Diagnostic (Explorers)
Endodontic Instruments


10



Page 86-87

Excavators

11



Page 88-89

Handles For Mouth Mirrors
Mouth Mirrors

12



Page 90-93

Spatulas



13

Needle Holders &
Saliva Ejectors


Page 94-100



14

Syringes & Scissors


Page 101-109



15

Dental Tweezers

Page 110-114



16

Rubber Dam Clamps


Page 115-117



17

Matrix Retainers &
Matrix Bands

Page 118-119








18

Crown Removers, Scissors,
Forceps, Amalgam Guns & Carriers

Page 120-122

TABLE OF CONTENTS

Page 123	Copper Ring Removing Pliers & Paper Articulating Forceps	19	
Page 124-127	Rongeurs	20	
Page 128-129	Bone Elevators, Files & Chisels	21	
Page 130-151	Pliers For Orthodontics, Mouth Gags, Tongue Depressors Lip and Cheek Retractors and Implantology Instrumetns	22	
Page 152-155	Impression Trays	23	
Page 156-157	Articulators	24	
Page 158-160	Hollowares	25	

Page
1-31

- Extracting Forceps - English Pattern







TK-1-5101
Fig 8
lower premolars



TK-1-5102
Fig 13
lower premolars



TK-1-5103
Fig 17
upper molars, right




TK-1-5104
Fig 18
upper molars, left




TK-1-5105
Fig 18 A
upper molars



TK-1-5106
Fig 19
upper wisdoms



TK-1-5107
Fig 20
lower wisdoms



TK-1-5108
Fig 21
lower molars



TK-1-5109
Fig 22
lower molars



TK-1-5110
Fig 22 A
Albrecht
lower molars



TK-1-5111
Fig 33 A
lower roots



TK-1-5112
Fig 33 B
lower roots



TK-1-5113
Fig 33 L
lower roots




TK-1-5114
Fig 33 L
lower roots




TK-1-5115
Fig 33 S
lower roots, for children




TK-1-5116
Fig 44
upper roots




TK-1-5117
Fig 45
lower roots



TK-1-5118
Fig 12
upper roots



TK-1-5119
Fig 51 A
upper roots




TK-1-5120
Fig 51 L
upper roots



Relax
TK-1-5121
Fig 51
upper roots




Relax
TK-1-5122
Fig 51 A
upper roots



Relax
TK-1-5123
Fig 51 L
upper roots



Relax
TK-1-5124
Fig 67
upper wisdoms



Relax
TK-1-5125
Fig 67 A
upper wisdoms




TK-1-5126
Fig 37
Upper incisors



TK-1-5127
Fig 39
Upper molars



TK-1-5128
Fig 39 R
Upper molars, right



TK-1-5129
Fig 39 L
Upper molars, left



TK-1-5130
Fig 515
Upper roots




Klein
TK-1-5131
Fig 6
lower molars



Klein
TK-1-5132
Fig 137
upper incisors



Klein
TK-1-5133
Fig 5
lower incisors




Klein
TK-1-5134
Fig 139
upper premolars

Extracting Forceps - English Pattern




TK-1-5135
Fig 43 A
Fractured or decayed
upper teeth



TK-1-5136
Fig 44
Upper roots



TK-1-5137
Fig 45
Lower roots



TK-1-5138
Fig 46
Lower roots




TK-1-5139
Fig 46 L
Lower roots




TK-1-5140
Fig 48
Upper roots



TK-1-5141
Fig 51
Upper roots



TK-1-5142
Fig 49
Upper roots



TK-1-5143
Fig 51 L
Upper roots



TK-1-5144
Fig 51 A
Upper roots

Extracting Forceps - English Pattern



TK-1-5145
Fig 51 LX
Upper roots



1046
Fig 52
Upper roots and
canines




TK-1-5147
Fig 54
For separating
upper molars



TK-1-5148
Fig 55
For separating
upper molars




TK-1-5149
Fig 56
For separating
lower molars



Extracting Forceps - English Pattern



TK-1-5155
Fig 66 L
Upper molars, left



TK-1-5156
Fig 66 R
Upper molars, right



TK-1-5157
Fig 67
Upper third molars



TK-1-5158
Fig 67 A
Upper third molars



TK-1-5159
Fig 68
Lower roots



TK-1-5160
Fig 67½ L
Upper third molars, left



TK-1-5161
Fig 67½ R
Upper third molars, right



TK-1-5162
Fig 73
Lower molars



TK-1-5163
Fig 74
Lower roots



TK-1-5164
Fig 74 M
Lower roots

Extracting Forceps - English Pattern




TK-1-515
Fig 74 N
Lower roots




TK-1-516
6 Fig 75
Lower canines and
premolars



TK-1-5167
Fig 79
Lower third molars



TK-1-5168
Fig 79 A
Lower third molars



TK-1-5169
Fig 79 C
Lower third molars



TK-1-5170
Fig 86 W
Lower third molars



TK-1-5171
Fig 86 1/2 L
Lower third molars, left



TK-1-5172
Fig 86 1/2 R
Lower third molars, right



TK-1-5173
Fig 86 C
Lower molars

Extracting Forceps - English Pattern



TK-1-5174
Fig 87
Lower decayed molars
or broken-down crowns




WITH SERRATED
TIPS

TK-1-5175
Fig 89
Upper molars, right




WITH SERRATED
TIPS

TK-1-5176
Fig 90
Upper molars, left



WITH SERRATED
TIPS

TK-1-5177
Fig 99
Broken-down lower molars




WITH SERRATED
TIPS

TK-1-5178
Fig 107
Upper canines



TK-1-5179
Fig 121
Upper third molars



WITH SERRATED
TIPS

TK-1-5180
Fig 145
Lower roots



WITH SERRATED
TIPS

TK-1-5181
Fig 151½
Upper roots



WITH SERRATED
TIPS

TK-1-5182
Fig 166
Universal extracting forceps for
upper teeth, in particular molars
and premolars



WITH SERRATED
TIPS

TK-1-5183
Fig 167
Universal extracting forceps for
lower teeth, in particular molars
and premolars

Extracting Forceps - English Pattern



TK-1-5184
Fig 168
Upper canines



TK-1-5185
Fig 169
Lower canines



TK-1-5186
Fig 346
Lower roots



TK-1-5187
Fig 349
Upper roots



TK-1-5188
Fig 359
Lower roots



TK-1-5189
Fig 1
Universal extracting forcep for
upper incisors and canines



TK-1-5190
Fig 2
Upper molars



TK-1-5191
Fig 3
Lower incisors, premolar
and roots




TK-1-5192
Fig 4
Lower molars




Extracting Forceps With Anatomically Shaped Handle



TK-1-5193
Fig 1
Upper incisors
and canines



TK-1-5194
Fig 2
Upper incisors
and canines



TK-1-5195
Fig 4
Lower incisors
and canines




TK-1-5196
Fig 7
Upper premolars



TK-1-5197
Fig 13
Lower premolars




TK-1-5198
Fig 17
Upper molars, right



WITH SERRATED TIPS

TK-1-5199
Fig 18
Upper molars, left



WITH SERRATED TIPS

TK-1-5200
Fig 21
Lower molars



WITH SERRATED TIPS

TK-1-5201
Fig 22
Lower molars



WITH SERRATED TIPS

TK-1-5202
Fig 33
Lower roots

Extracting Forceps With Anatomically Shaped Handle



TK-1-5203
Fig 45
Lower roots



WITH SERRATED
TIPS



TK-1-5204
Fig 51 A
Upper roots

WITH SERRATED
TIPS



TK-1-5205
Fig 67 A
Upper third
molars

WITH SERRATED
TIPS



TK-1-5206
Fig 74
Lower roots

WITH SERRATED
TIPS



TK-1-5207
Fig 79
Lower third
molars

WITH SERRATED
TIPS



TK-1-5208
Fig 86 C
Lower molars

WITH SERRATED
TIPS



TK-1-5209
Fig 359
Lower roots

WITH SERRATED
TIPS

Extracting Forceps With Anatomically Shyped Handle




TK-1-5210
Fig 89
Upper molars, right



WITH SERRATED
TIPS

TK-1-5211
Fig 90
Upper molars, left




WITH SERRATED
TIPS

TK-1-5212
Fig 145
Lower roots



WITH SERRATED
TIPS

TK-1-5213
Fig 151 1/2
Upper roots



WITH SERRATED
TIPS

TK-1-5214
Fig 349
Upper roots



TK-1-5215
Fig 5
Lower third molars



TK-1-5216
Fig 6
Separating forcep



TK-1-5217
Fig 105
Upper molars



TK-1-5218
Fig 16
Lower molars



TK-1-5219
Fig 17
Lower molars



TK-1-5220
Fig 18 L
Upper molars, left



TK-1-5221
Fig 18 R
Upper molars, right



TK-1-5222
Fig 23
Lower molars



TK-1-5223
Fig 53 R
Upper molars, right



TK-1-5224
Fig 53 L
Upper molars, left

Extracting Forceps - American Pattern



TK-1-5225
Fig 65
Upper incisors and roots



TK-1-5226
Fig 69
Upper and lower roots



TK-1-5227
Fig 88 R
Upper molars, right



TK-1-5228
Fig 88 L
Upper molars, left



TK-1-5229
Fig 101
Deciduous teeth



TK-1-5230
Fig 150
Upper incisors, premolars, roots



TK-1-5231
Fig 150 S
Upper teeth (for children)



TK-1-5232
Fig 151
Upper incisors, premolars, roots



TK-1-5233
Fig 151S
Lower teeth (for children)



TK-1-5234
Fig 222
Lower third molars




TK-1-5235
Fig 1
upper incisors and bicuspid



TK-1-5236
Fig 1 B
upper roots, anterior



Physick
TK-1-5237
Fig 5
lower wisdoms



TK-1-5238
Fig 6
lower molars, splitting forceps



TK-1-5239
Fig 10 H
upper molars

Extracting Forceps For Children - English Pattern



TK-1-5240
Fig 135
Lower premolars

WITH SERRATED
TIPS



TK-1-5241
Fig 225
Lower molars



TK-1-5242
Fig 305
Upper incisors
and roots



TK-1-5243
Fig 33 S
Lower roots

WITH SERRATED
TIPS



TK-1-5244
Fig 38
Lower incisors
and canines

WITH SERRATED
TIPS

Extracting Forceps For Children - Klein Pattern



TK-1-5245
Fig 1
Upper incisors



TK-1-5246
Fig 2
Upper premolars



TK-1-5247
Fig 3
Upper molars



TK-1-5248
Fig 4
Upper roots



TK-1-5249
Fig 5
Lower incisors



TK-1-5250
Fig 6
Lower molars



TK-1-5251
Fig 7
Lower roots



TK-1-5252
Fig 1
Lower roots



TK-1-5253
Fig 2
Upper roots



TK-1-5255
Fig 2
Roots fragments



TK-1-5254
Fig 1
Roots fragments


Page 32–50

- Root Elevators



lower jaw, left side

lower jaw, right side



TK-2-5101



TK-2-5102

Root elevator, spiral-shaped, for extremely facilitated extraction of the molars in the lower jaw. The slender, conical points allow deep penetration between tooth and gum. Owing to the hook-shaped working parts, slipping of the elevator is avoided.

Root Elevators



Bein



Bein
TK-2-5103
Fig 1



Bein
TK-2-5104
Fig 2



Bein
TK-2-5105
Fig 0



Bein
TK-2-5106
Fig 3



Bein
TK-2-5107
Fig 4



Bein
TK-2-5108
Fig 5



Bein
TK-2-5109
Fig 6



Seldin



**Seldin
TK-2-5110
Fig 1L**



**Seldin
TK-2-5111
Fig 1R**



**Seldin
TK-2-5112
Fig 2**



**Seldin
TK-2-5113
Fig 3**



**Seldin
TK-2-5114
Fig 4L**



**Seldin
TK-2-5115
Fig 4R**



**Cryer
TK-2-5116
Fig 3**



**Cryer
TK-2-5117
Fig 4**



**Cryer
TK-2-5118
Fig 12**



**Cryer
TK-2-5119
Fig 27**



**Cryer
TK-2-5120
Fig 28**



**Cryer
TK-2-5121
Fig 27 S**



**Cryer
TK-2-5122
Fig 28 S**



Apical



Apical
TK-2-5123
Fig 1



Apical
TK-2-5124
Fig 2



Apical
TK-2-5125
Fig 3



Apical
TK-2-5126
Fig 4



Cryer



**Cryer
TK-2-5127
Fig 1**



**Cryer
TK-2-5128
Fig 2**



Coupland



**Coupland
TK-2-5129
Fig 1**



**Coupland
TK-2-5130
Fig 2**




**Coupland
TK-2-5131
Fig 3**



TK-2-5132
Fig 1



TK-2-5133
Fig 2



TK-2-5134
Fig 3



Flohr
TK-2-5135
Fig 1



Flohr
TK-2-5136
Fig 2



Flohr
TK-2-5137
Fig 3

Root Elevators



Kopp
TK-2-5138
Fig 1



Kopp
TK-2-5139
Fig 2



Kopp
TK-2-5144
Fig 1



Kopp
TK-2-5145
Fig 2



Kopp
TK-2-5140
Fig 3



Kopp
TK-2-5141
Fig 4



Kopp



Kopp
TK-2-5146
Fig 3




Kopp
TK-2-5147
Fig 4



Kopp
TK-2-5142
Fig 5



Kopp
TK-2-5143
Fig 6



Hylin
TK-2-5148
Fig 1




Hylin
TK-2-5149
Fig 2




Bein



Bein
TK-2-5150
Fig 1



Bein
TK-2-5151
Fig 2



Bein
TK-2-5152
Fig 3



Bein
TK-2-5153
Fig 4




**Warwick-James
TK-2-5154
Fig 1**



**Warwick-James
TK-2-5155
Fig 2**



**Warwick-James
TK-2-5156
Fig 3**



**London
Hospital
TK-2-5157
Fig 1**




**London
Hospital
TK-2-5158
Fig 2**



**London
Hospital
TK-2-5159
Fig 3**



**Ott
TK-2-5160**



Pott



**Pott
TK-2-5162
Fig 1**



**Pott
TK-2-5163
Fig 2**



**TK-2-5164
Fig 1**



**TK-2-5165
Fig 2**



Mead
TK-2-5166
MD 1



Mead
TK-2-5167
MD-5 A



Mead
TK-2-5168
MD-6 A



Mead
TK-2-5169
MD-7 A



Mead
TK-2-5170
MD-8



Lecluse
TK-2-5171



**Chompret
TK-2-5172
Fig 1**



**Chompret
TK-2-5173
Fig 2**



**Chompret
TK-2-5174
Fig 3**



**Chompret
TK-2-5175
Fig 4**



**Apical
TK-2-5176**

Root Elevators Root Tip Elevators



TK-2-5177
Fig 1



TK-2-5178
Fig 2



TK-2-5179
Fig 3



Heidbrink
TK-2-5180
Fig 1



Heidbrink
TK-2-5181
Fig 2



Heidbrink
TK-2-5182
Fig 3



Winter
TK-2-5183
Fig 1L



Winter
TK-2-5184
Fig 1R



Winter
TK-2-5185
Fig 2L



Winter
TK-2-5186
Fig 2R



Winter



Winter
TK-2-5187
Fig 11 L



Winter
TK-2-5188
Fig 11 R



Winter
TK-2-5189
Fig 12 L



Winter
TK-2-5190
Fig 12 R



Winter
TK-2-5191
Fig 14 L



Winter
TK-2-5192
Fig 14 R



Barry
TK-2-5193
Fig 1



Barry
TK-2-5194
Fig 2



Barry
TK-2-5195
Fig 3



Barry
TK-2-5196
Fig 4



Barry



Barry
TK-2-5197
Fig 11



Barry
TK-2-5198
Fig 12



Barry
TK-2-5199
Fig 17



Barry
TK-2-5200
Fig 18



Cryer
TK-2-5201
Fig 26



Seldin
TK-2-5202
Fig 27



Cryer
TK-2-5203
Fig 28



Cryer
TK-2-5204
Fig 29



Heidbrink
TK-2-5205
Fig 22



Heidbrink
TK-2-5206
Fig 23



TK-2-5207



Kopp
TK-2-5208
Fig 31



Kopp
TK-2-5209
Fig 32



Kopp
TK-2-5210
Fig 33



Kopp
TK-2-5211
Fig 34



Kopp
TK-2-5212
Fig 35



Kopp
TK-2-5213
Fig 44



TK-2-5214



Heidbrink
TK-2-5215
Fig 10



Winter
TK-2-5216
Fig 11



Winter
TK-2-5217
Fig 12



Barry
TK-2-5218
Fig 17



Barry
TK-2-5219
Fig 18



TK-2-5220
Fig 20



TK-2-5221

Root Elevators according to Dr. Sader,
with interchangeable working tips



TK-2-5222



TK-2-5223




TK-2-5224

Page 51–55

- Filling Instruments





TK-2-5226 Williams TK-2-5227 Williams TK-2-5228 Williams TK-2-5229 Williams


TK-2-5234 White TK-2-5235 White TK-2-5236 White


TK-2-5225 Williams

TK-2-5230 Williams TK-2-5231 Williams TK-2-5232 Williams

TK-2-5233 White

TK-2-5237 White TK-2-5238 White TK-2-5239 White







TK-2-5268
Heidemann
Fig 1



TK-2-5269
Heidemann
Fig 2



TK-2-5270
Heidemann
Fig 3



TK-2-5271



TK-2-5272
Fig 146




TK-2-5274
Fig 148




TK-2-5273
Fig 147




TK-2-5275
Fig 149




TK-2-5276
Fig 158



TK-2-5277
Fig 159



TK-2-5278
Woodson
Fig 160



Page 56–59

- Dental Probes



With Point Set In



TK-2-5295



TK-2-5296



TK-2-5297



TK-2-5298



TK-2-5299



TK-2-5300



TK-2-5301



TK-2-5302



TK-2-5303



TK-2-5304



TK-2-5305




TK-2-5306



TK-2-5307
Woodson

Polishing Instruments





- New marking process.
- Laser system marking.
- No more paint peeling off.
- Permanent marking.
- Antiglare, low depth marking.
- Maximum tip resistance against breaking.



TK-2-5322
Fig CP-C



TK-2-5323
Fig CP-11



TK-2-5324
Fig CP-12



TK-2-5325
Fig CP-12S
O.M.S



TK-2-5326
Williams



TK-2-5327
Merrit A



TK-2-5328
Merrit B



TK-2-5329
Goldman
Fox



TK-2-5330
University of
Michigan



TK-2-5331
26 G

Page 60–67

- Periodontia Instruments, Mallets



Files

Chisels



TK-3-5101
Sugarman
Fig. 1S/2S



TK-3-5102
Sugarman
Fig. 3S/4S



TK-3-5103
Schluger
Fig. 9/10



TK-3-5104
Buck
Fig. 11/12



TK-3-5105
Ochsenein
Fig. 3




TK-3-5106
Ochsenein
Fig. 4



TK-3-5107
Rhodes
Back - Action
Fig. 36/37

Files



TK-3-5108
Sugarman
Fig. 1S/2S

TK-3-5109
Sugarman
Fig. 3S/4S

TK-3-5110
Schluger
Fig. 9/10



TK-3-5111
Buck
Fig. 11/12



TK-3-5112
Ochsenbein
Fig. 3



TK-3-5113
Ochsenbein
Fig. 4



TK-3-5114
Rhodes
Back - Action
Fig. 36/37

Gingivectomy knives



TK-3-5115
Goldman-Fox
Fig. 7



TK-3-5116
Goldman-Fox
Fig. 8



TK-3-5117
Goldman-Fox
Fig. 9



TK-3-5118
Goldman-Fox
Fig. 11



TK-3-5119
Buck
Fig. 3/4



TK-3-5120
Buck
Fig. 5/6

Curettes and Scalers



TK-3-5121
Curette
Fig. HU 1/2



TK-3-5122
Columbia
Fig. 2R/2L



TK-3-5123
Columbia
Fig. 4R/4L



TK-3-5124
Columbia
Fig. 13/14



TK-3-5125
Fig. 411/412



TK-3-5126
Fig. 413/414



TK-3-5127
Zerfing-chisel

Gingivectomy knives



TK-3-5128
Goldman - Fox
Fig. 7



TK-3-5129
Goldman - Fox
Fig. 8



TK-3-5130
Goldman - Fox
Fig. 9



TK-3-5131
Goldman - Fox
Fig. 11



TK-3-5132
Buck
Fig. 3/4



TK-3-5133
Buck
Fig. 5/6

Curettes and Scalers



TK-3-5134
Curette
Fig. HU 1/2



TK-3-5135
Columbia
Fig. 2R/2L



TK-3-5136
Columbia
Fig. 4R/4L



TK-3-5137
Columbia
Fig. 13/14



TK-3-5138
Fig. 413/412




TK-3-5139
Fig. 411/414



TK-3-5140
Zerfing-chisel

Bone Chisels and Gouges



TK-3-5141
Ochsenbein
Fig. 1

TK-3-5142
Ochsenbein
Fig. 2

TK-3-5143
Gardner
Fig. 1

TK-3-5144
Gardner
Fig. 2

TK-3-5145
Gardner
Fig. 5

TK-3-5146
Gardner
Fig. 5R


TK-3-5147
Gardner
Fig. 6

TK-3-5148
Wakefield
Fig. 1W


TK-3-5149
Wakefield
Fig. 2W

TK-3-5150
Buckley
Fig. 3

Bone Chisels and Gouges




Martelli chirurgici Mallet



Page 68–74

▪ Currettes Scalers





TK-3-5170
Gracey
Fig. 1/2



TK-3-5171
Gracey
Fig. 3/4



TK-3-5172
Gracey
Fig. 5/6



TK-3-5173
Gracey
Fig. 7/8



TK-3-5174
Gracey
Fig. 9/10



TK-3-5175
Gracey
Fig. 11/12



TK-3-5176
Gracey
Fig. 13/14




TK-3-5177
Gracey
Mesial Distal
Fig. 11/14



TK-3-5178
Gracey
Mesial Distal
Fig. 12/13

Gracey Curettes AF Type



TK-3-5179
Gracey
AF/type
Fig. 1/2



TK-3-5180
Gracey
AF/type
Fig. 3/4



TK-3-5181
Gracey
AF/type
Fig. 5/6



TK-3-5182
Gracey
AF/type
Fig. 7/8




TK-3-5183
Gracey
AF/type
Fig. 11/12




TK-3-5184
Gracey
AF/type
Fig. 13/14



TK-3-5185
Gracey
AF/type
Fig. 15/16




TK-3-5187 Gracey AF/type Fig. 1/2	TK-3-5188 Gracey AF/type Fig. 3/4	TK-3-5189 Gracey AF/type Fig. 5/6	TK-3-5190 Gracey AF/type Fig. 7/8	TK-3-5191 Gracey AF/type Fig. 11/12	TK-3-5192 Gracey AF/type Fig. 13/14	TK-3-5193 Gracey AF/type Fig. 15/16
---	---	---	---	---	---	---



TK-3-5186 Gracey AF/type	TK-3-5194 Gracey MF/type Fig. 1/2	TK-3-5195 Gracey MF/type Fig. 3/4	TK-3-5196 Gracey MF/type Fig. 5/6	TK-3-5197 Gracey MF/type Fig. 7/8	TK-3-5198 Gracey MF/type Fig. 11/12	TK-3-5199 Gracey MF/type Fig. 13/14	TK-3-5200 Gracey MF/type Fig. 15/16
---------------------------------------	---	---	---	---	---	---	---

Gracey Curettes AF Type



TK-3-5201
Gracey
MF/type
Fig. 1/2



TK-3-5202
Gracey
MF/type
Fig. 3/4



TK-3-5203
Gracey
MF/type
Fig. 5/6



TK-3-5204
Gracey
MF/type
Fig. 7/8



TK-3-5205
Gracey
MF/type
Fig. 11/12



TK-3-5206
Gracey
MF/type
Fig. 13/14



TK-3-5207
Gracey
MF/type
Fig. 15/16



TK-3-5208
Goldman - Fax
Fig. 1



TK-3-5209
Goldman - Fax
Fig. 2




TK-3-5210
Goldman - Fax
Fig. 3



TK-3-5211
Goldman - Fax
Fig. 4



TK-3-5212
Goldman - Fax
Fig. 21




TK-3-5213
Goldman - Fax
Fig. 5



TK-3-5214
Mc Call
Fig. 13S/14S



TK-3-5215
Younger-Good
Fig. 7/8





TK-3-5228
Hemingway
Fig. 1



TK-3-5229
Hemingway
Fig. 2



TK-3-5230
Hemingway
Fig. 3



TK-3-5231
Lucas
Fig. 85



TK-3-5232
Lucas
Fig. 86




TK-3-5233
Lucas
Fig. 87


Page 75–77

- Modelling Instruments & Measuring Instruments





TK-4-5101 Le Cron TK-4-5102 Le Cron TK-4-5103 Le Cron TK-4-5104 MSY TK-4-5105 MSY TK-4-5106 Le Cron TK-4-5107 Hylin TK-4-5108 Hylin TK-4-5109 Hylin TK-4-5110 Zahle TK-4-5111 Zahle TK-4-5112 Zahle



TK-4-5113 Gritman TK-4-5114 Beale TK-4-5115 Roach TK-4-5116 Vehe TK-4-5117 P.K. Thomas Fig. 1 TK-4-5118 P.K. Thomas Fig. 2 TK-4-5119 P.K. Thomas Fig. 3 TK-4-5120 P.K. Thomas Fig. 4 TK-4-5121 P.K. Thomas Fig. 5



TK-4-5122
Boley




TK-4-5123
Zuricher
130 mm



TK-4-5124
Vernier
12 cm



TK-4-5125
Vernier




TK-4-5126
Carroll



TK-4-5127
Castroviejo's
Scleral Marker



TK-4-5128
Iwanson



TK-4-5129
Albee



TK-4-5130
18 cm
7 1/8"



TK-4-5131
18 cm
7 1/8"



TK-4-5132
17 cm
6 3/4"
Gritman



TK-4-5133
17 cm
6 3/4"
Medicon



TK-4-5134
22 cm
8 3/4"



TK-4-5135
Antkoviak



TK-4-5136
Ermert



TK-4-5137
Ermert

Page 78–79

▪ Wax Knives/plaster





TK-4-5138
Fahnenstock
Fig. 1
12,5cm



TK-4-5139
Fahnenstock
Fig. 2
17cm



TK-4-5140
Lessmann
Fig. 1
12,5cm



TK-4-5141
Lessmann
Fig. 2
17cm



TK-4-5142
Con scodellino
Fig. 1
12,5cm



TK-4-5143
Con scodellino
Fig. 2
17cm



TK-4-5144
Fig. 1
17 cm



TK-4-5145
Gritman
Fig. 2
17 cm



TK-4-5146
Fig. 3
16 cm



TK-4-5147
Fig. 4
16 cm

Page 80–85

▪ Diagnostic, Explorer, Endodontic Instruments



Explorers



TK-4-5148
Fig 23



TK-4-5149
Fig 6



TK-4-5150
Fig 8



TK-4-5151
Fig 9



TK-4-5152
Fig 17




TK-4-5153
Fig 3
Rein



TK-4-5154
Fig 23/17



TK-4-5155
Fig 17



TK-4-5156
Fig.1
Spreader

TK-4-5157
Fig.2
Spreader

TK-4-5158
Fig.3
Spreader

TK-4-5159
Fig.4
Spreader


TK-4-5110
Fig.5
Spreader

TK-4-5111
Fig.6
Spreader

TK-4-5112
Fig.7
Spreader


Double-Ended Root Canal Explorers Plugger/Spreaders

Ø 0.4 Ø 0.4 Ø 0.6



Spreaders


Ø 0.6 Ø 0.5



Periodontal Pocket Probes



Root-canal pluggers




Endodontic set




TK-4-5132
MeasuringPlate

Root-canal pluggers

Optimum Flexibility
Very Fine Points



Root-canal explorers



TK-4-5143
Fig.6/1

TK-4-5144
Fig.6/17

TK-4-5145
Fig.6/23

Root-canal excavators



TK-4-5146
Fig.31L

TK-4-5147
Fig.32L

TK-4-5148
Fig.33L

Page 86–87

- Excavators



1 mm



TK-4-5149
Fig 5

1.2 mm



TK-4-5150
Fig 14

1 mm



TK-4-5151
Fig 17

1.2 mm



TK-4-5152
Fig 18

1.6 mm



TK-4-5153
Fig 19

2 mm



TK-4-5154
Fig 20

1 mm



TK-4-5155
ASH
Fig.243

1,4 mm



TK-4-5156
ASH
Fig.G2

0,8 mm



TK-4-5157
Fig.1

0,8 mm



TK-4-5158
Fig.2

1 mm



TK-4-5159
Fig.11

1 mm



TK-4-5160
Fig.12

1,2 mm



TK-4-5161
Fig.21

1,2 mm



TK-4-5162
Fig.22

1,6 mm



TK-4-5163
Fig.31

1,6 mm



TK-4-5164
Fig.32

2 mm



TK-4-5165
Fig.41

2 mm



TK-4-5166
Fig.42

Page 88–89

- Handles For Mouth Mirrors



Handles for Mouth Mirrors and Mouth Mirrors

Set of 5 pieces

TK-5-5101 TK-5-5102 TK-5-5103

Hollow, stainless steel

Plane

TK-5-5104
Fig.6
Ø 26 mm

TK-5-5105
Fig.5
Ø 24 mm

TK-5-5106
Fig.4
Ø 22 mm

TK-5-5107
Fig.4
Rhodium
Ø 22 mm

Front surface mouth mirrors



Magnifying

TK-5-5108
Fig.6
Ø 26 mm

TK-5-5109
Fig.5
Ø 24 mm

TK-5-5110
Fig.4
Ø 22 mm

TK-5-5111
Fig.5
Rhodium
Ø 24 mm



Page 90–93

- Spatulas





TK-5-5126




TK-5-5127



TK-5-5128



TK-5-5129



TK-5-5130




TK-5-5131



TK-5-5132




TK-5-5133





Spatulas for Plaster and Alginate



Page 94–100

- Needle Holders





TK-5-5150
Mayo-Hegar
16 cm



TK-5-5151
Mayo-Hegar TC
16 cm



TK-5-5152
Mayo-Hegar
18 cm



TK-5-5153
Mayo-Hegar TC
18 cm



TK-5-5154
Mini-Ryder
18 cm



TK-5-5155
Mini-Ryder
15 cm



TK-5-5156
Derf
12 cm



TK-5-5157 Hegar-
Baumgarther 14 cm



TK-5-5158
Crile-Wood
15 cm



TK-5-5159
Crile-Wood TC
15 cm



TK-5-5160
Adson
17.5 cm



TK-5-5161
Olsen-Hegar
14 cm



TK-5-5162
Baumgarther
14.5 cm



TK-5-5163
Gillies
16 cm



TK-5-5164
Gillies
16 cm



TK-5-5165 Hegar-
Baumgarther 14 cm



TK-5-5166
Toennis
18 cm



TK-5-5167
Toennis
18 cm



TK-5-5168
Lichtenberg
17 cm



TK-5-5169
Lichtenberg
17 cm



TK-5-5170
Eiselberg-Mathieu
19.5 cm



TK-5-5171
Boynton
12 cm



TK-5-5172
Arruga
16 cm



TK-5-5173
Boynton
12 cm



TK-5-5174
Arruga
16 cm



TK-5-5175
Castroviejo
14 cm



TK-5-5176
Castroviejo
14 cm



TK-5-5177
Castroviejo
14/60 cm



TK-5-5178
Castroviejo
18 cm



TK-5-5179
Castroviejo
18 cm



TK-5-5180
Castroviejo
14 cm



TK-5-5181
Castroviejo
14 cm



TK-5-5182
Castroviejo
14/60 cm



TK-5-5183
Castroviejo
18 cm



TK-5-5184
Castroviejo
18 cm



TK-5-5185
Mathieu
14 cm



TK-5-5186
Mathieu
14 cm



TK-5-5187
Mathieu
14 cm



TK-5-5188
Mathieu
14 cm



TK-5-5189
Mathieu
17 cm



TK-5-5190
Mathieu
17 cm



TK-5-5191
Mathieu
14 cm



TK-5-5192
Micro-Barraquer
18 cm

Needle Holders & Stainless Steel Saliva Ejectors



The tips have tungsten carbide inserts for better handling.



TK-5-5193
Mathieu
17 cm



TK-5-5194
Mathieu
17 cm



TK-5-5195
Micro-Barraquer
14 cm



TK-5-5196
Castroviejo
13 cm



TK-5-5197
Castroviejo-Barraquer
13 cm



TK-5-5198
Castroviejo
13 cm



TK-5-5199
Castroviejo
13 cm



TK-5-5200
18.03mm,cm






TK-5-5201
18.04mm,cm



TK-5-5202
18.05mm,cm



TK-5-5206
Fig.6FR
Fraizer
15.02mm,cm



Page 101–109

- Syringes and Scissors





TK-5-5222
Water Syringe complete with
1 cannula (Lure-Lock)



TK-5-5223
Water Syringe complete with
1 cannula (Lure-Lock)



TK-5-5224
Metal water syringe, self-filling,
with 1 cannula (Lure-Lock)



TK-5-5225
Sparrel Barrel with
plunger suitable for
syringes



TK-5-5226 Cannula,
only, with
Luer-Lock connector
suitable for syringes



TK-5-5227
Aspirating Syringe



TK-5-5228
Position Gauge



TK-5-5229
Breech Loading Syringe



TK-5-5230
Aspirating & Breech Syringe



TK-5-5231
Band Adaptor



TK-5-5232
Endo Finger Ruler
With Adjustable Slide



TK-5-5233
Aspirating Syringe



TK-5-5234



TK-5-5236



TK-5-5235




TK-5-5239
Tip Only, short



TK-5-5240
Tip Only, long



TK-5-5241



TK-5-5237 1,8 ml

TK-5-5238 2,2 ml

Cartridge Syringe complete
With 2 tips and without needles



TK-5-5243



TK-5-5242



TK-5-5244
Cork screw tip



TK-5-5245
Cork screw tip



TK-5-5246
Folding type, cork
screw tip



TK-5-5247
Folding type, cork
screw tip



TK-5-5248



TK-5-5249
Cork screw tip



TK-5-5250
Spare cap with
metric thread



TK-5-5251
With Whitworth
thread



AS-5252
Autoclavable Plastic Protection sleeve



TK-5-5253
Brass syringe with stainless steel moving parts









TK-5-5254
Brass syringe with stainless
steel moving parts

Operating Scissors




Suture Scissors | Tijeras para Hilos | Fadenscheren





Scissors, One Blade Serrated


- One blade is with fine serrations to avoid tissue or ligature slippage.
- Manufactured from high quality stainless steel to provide longer-lasting cutting edge.
- Single black ring for easy identification.




Iris
11.5 cm




TK-5-5301





TK-5-5302



TK-5-5303




TK-5-5304




Goldman-Fox
13 cm

TK-5-5305




TK-5-5306
Goldman-Fox
13 cm



TK-5-5307



TK-5-5308



Kelly
16 cm



TK-5-5309
Locklin
16 cm



TK-5-5310
Locklin
16 cm

Super-Cut

- Cutting Blade horned to "razor" sharpness with opposing blade finished with fine serrations to avoid tissue slippage.
- Extremely smooth cutting with less trauma.
- Manufactured from high quality stainless steel to provide longer-lasting cutting edge.
- Single gold plated ring for easy identification.




TK-5-5311



TK-5-5312


Iris
11.5 cm



TK-5-5313
La Grange
11.5 cm



TK-5-5314




TK-5-5315


Goldman-Fox
13 cm




Iris
11.5 cm




Metzenbaum
14.5 cm




Kelly
16 cm



TK-5-5322
La Grange
11.5 cm




Goldman-Fox
13 cm



Strabismus
11.5 cm


Micro Scissors



TK-5-5327
15 cm

TK-5-5328
15 cm

Micro Scissor
Stainless Steel



TK-5-5329
18 cm

TK-5-5330
18 cm

Micro Scissor
Stainless Steel




TK-5-5331
18 cm

Micro Scissor
Stainless Steel
for ligature, one blade serrated

Micro Scissors with Titanium Stainless Steel Tips




TK-5-5332
18 cm



TK-5-5333
18 cm

TK-5-5334
18 cm



TK-5-5335
15.5 cm

TK-5-5336
15.5 cm



Page 110–114

- Dental Tweezers





smooth

TK-5-5337
College
15 cm



serrated

TK-5-5338
College
15 cm



smooth

TK-5-5339
College
15 cm



serrated

TK-5-5340
College
15 cm



smooth

TK-5-5341
College
15 cm



serrated

TK-5-5342
College
15 cm



smooth

TK-5-5343
College
15 cm



serrated

TK-5-5344
College
15 cm





TK-5-5345
Meriam
16 cm



smooth



TK-5-5346
Meriam
16 cm



serrated



TK-5-5347
Meriam
16 cm



smooth



TK-5-5348
Meriam
16 cm



serrated



TK-5-5349
Meriam
16 cm



smooth



TK-5-5350
Meriam
16 cm



serrated



TK-5-5351
Meriam
16 cm




smooth



TK-5-5352
Meriam
16 cm



serrated




TK-5-5353
College
14.5 cm, 5³/₄"

TK-5-5354
Flagg
16 cm, 6¹/₄"

TK-5-5355
Flagg
16 cm, 6¹/₄"

TK-5-5356
Landgraf
15 cm, 6"

TK-5-5357
15,5 cm, 6¹/₈"



TK-5-5358
15.5 cm, 6¹/₈"

TK-5-5359
London-College
15 cm, 6"

TK-5-5360
London-College
15 cm, 6"

TK-5-5361
London-College
15 cm, 6"

TK-5-5362
London-College
15 cm, 6"



TK-5-5363
CK1
15 cm, 6"

TK-5-5364
CK2
15 cm, 6"

TK-5-5365
GF2
16 cm, 6 1/4"

TK-5-5366
Gf1
16 cm, 6 1/4"

TK-5-5367
Kopp
15.5 cm, 6 1/8"



TK-5-5368
Kopp
15.5 cm, 6 1/8"

TK-5-5369
15 cm, 6"


TK-5-5370
16.5 cm, 6 1/2"

TK-5-5371
15.5 cm, 6 1/8"


TK-5-5372
15.5 cm, 6 1/8"



TK-5-5373
London-College
15 cm, 6"




TK-5-5374
London-College
15 cm, 6"




TK-5-5375
Meriam
16 cm, 6 1/4"



TK-5-5376
Allen
18 cm, 7 1/8"



TK-5-5377
Albrecht
16 cm, 6 1/4"



TK-5-5378
Albrecht
16 cm, 6 1/4"



TK-5-5379
16 cm, 6 1/4"



TK-5-5380
16 cm, 6 1/4"




TK-5-5381
16 cm, 6 1/4"

Page 115-117


- Rubber Dam Clamps



Rubber Dam Clamp Forceps



CLAMPS (PLIERS)





TK-5-5392
Winged
Labial Clamp



TK-5-5393 Winged
Edged Labial Clamp



TK-5-5394
Winged
Cervical Clamp



TK-5-5395
Winged
Labial Clamp

Premolars



TK-5-5396
Winged



TK-5-5397
Winged



TK-5-5398
Wingless

For Small Teeth, Elongated or With Narrow Collars



TK-5-5399
Winged



TK-5-5400
Winged



TK-5-5401
Wingless

For General Use With Adjustable Jaws According
To Gingival Line Shapes



TK-5-5402
Winged



TK-5-5403
Edged, Winged



TK-5-5404
Wingless

For Premolars, With Elongated Collar



TK-5-5405
Winged



TK-5-5406
Edged, Winged



TK-5-5407
Wingless

For Premolars, With Elongated Collar



TK-5-5408
Wingless



TK-5-5409
Wingless

Rubber Dam Clamps



TK-5-5410
Wingless

TK-5-5411
Wingless

For Small Lower Molars



TK-5-5412
Winged

TK-5-5413
Wingless

For Small Lower Molars With Round clamps



TK-5-5414
Winged

TK-5-5415
Wingless

For Square Lower Molars



TK-5-5416
Winged

TK-5-5417
Wingless

For Round Upper Molars



TK-5-5418
Winged

TK-5-5419
Wingless

For Molar Roots



TK-5-5420
Winged

TK-5-5421
Wingless

For Partially Broken Molars



TK-5-5422
Winged

TK-5-5423
Edged Winged

TK-5-5424
Wingless

For Round Molars, Partially Broken, Undersized or Inclined



TK-5-5425
Winged

TK-5-5426
Wingless



TK-5-5427
Wingless

TK-5-5428
Wingless



TK-5-5429
Winged

TK-5-5430
Winged

Page 118–122

- Matrix Retainer
- Matrix Bands



Matrix Retainers and Matrix Bands



TK-6-5101
Ivory
63mm



TK-6-5102
Tofflemire-Universal



TK-6-5103
Tofflemire-Junior



TK-6-5104
Tofflemire-Senior



TK-6-5105
Ivory 8
5 mm



TK-6-5106
Ivory 8
6 mm



TK-6-5107
Ivory 8
7 mm



TK-6-5108
Ivory 8 N



TK-6-5109
Didner 8 P

1750-5 mm x 100 cm
1750-6 mm x 100 cm
1750-7 mm x 100 cm



TK-6-5110



TK-6-5111



TK-6-5112



TK-6-5113



TK-6-5114

For Molars and Premolars



TK-6-5115




TK-6-5116



TK-6-5117




7
TK-6-5118
Molar
Fig 7



8
TK-6-5119
Molar
Fig 8



9
TK-6-5120
Molar
Fig 9



10
TK-6-5121
Molar
Fig 10



1
TK-6-5122
Tofflemire



2
TK-6-5123
Tofflemire



3
TK-6-5124
Tofflemire



11
TK-6-5125
Bicuspid



12
TK-6-5126
Bicuspid



13
TK-6-5127
Bicuspid



1
TK-6-5128
Bicuspid



2
TK-6-5129
Bicuspid



3
TK-6-5130
Bicuspid



TK-6-5134



1
TK-6-5131
Molar



2
TK-6-5132
Molar







3
TK-6-5133
Molars


Page 120–122

- Crown Removers
- Scissors
- Forceps
- Amalgum Guns
- Carries









TK-6-5156 Medium-Large 20 cm
TK-6-5157 Medium-Jumbo 20 cm
TK-6-5158 Large-Jumbo 20 cm
TK-6-5159 Jumbo-Jumbo 20 cm
TK-6-5160 Medium 15.5 cm
TK-6-5161 Large 15.5 cm
TK-6-5162 Jumbo 15.5 cm



TK-6-5163 Straight 15 cm
TK-6-5164 15 cm
TK-6-5165 15 cm



TK-6-5166
Amalgam Well
Non Skid Base
Ø 3.5 mm

Page 123

- Copper Ring Removing Plierspaper Articulator Forceps



Copper Ring Removing Pliers, Paper Articulator Forceps

Copper ring removing pliers



TK-6-5167
Bohm
11,5 cm



TK-6-5168
Bohm
14,5 cm



TK-6-5169
Furrer
16 cm



TK-6-5170
Furrer
16 cm



TK-6-5171
Miller
15 cm



TK-6-5172
Miller
15 cm



TK-6-5173
Kelly
14 cm



TK-6-5174
Halstead-Mosquito
12 cm




TK-6-5175
Micro-Mosquito
10 cm


Page 124–127

▪ Rongeurs







TK-6-5176
Friedmann
14 cm, 5 1/2"



TK-6-5177
14 cm, 5 1/2"




TK-6-5178
Mead
14 cm, 5 1/2"



TK-6-5179
Mead
14 cm, 5 1/2"



TK-6-5180
Mead
17 cm, 6 3/4"



TK-6-5181
Schell
17 cm, 6 3/4"



TK-6-5182
Luer
15 cm



TK-6-5183
Friedman
14 cm



TK-6-5184
Mini Friedman
12 cm



TK-6-5185
Goldman-Fax
10 cm



TK-6-5186
Blumenthal
15 cm



TK-6-5187
Blumenthal 45
14,5 cm




TK-6-5188
Blumenthal 90
14 cm



TK-6-5189
Bohler
15 cm




TK-6-5190
Kerrison
15 cm



TK-6-5191
Blumenthal
16 cm, 6 1/4"



TK-6-5192
Blumenthal
16 cm, 6 1/4"



TK-6-5193
Blumenthal
16 cm, 6 1/4"



TK-6-5194
Kleinert-Kutz
14 cm, 5 1/2"




TK-6-5195
Ruskin-Mini
15 cm



TK-6-5196
Beyer
18 cm, 7 1/8"



TK-6-5197
Zaufal-Jansen
18 cm, 7 1/8"



TK-6-5198
Niro
No. 5 S, 14 cm, 5 1/2"



TK-6-5199
Cleveland
No. 5 S, 17 cm, 6 3/4"



TK-6-5200
No. 18, 15 cm, 6"



TK-6-5201
Cleveland
16 cm, 6 1/4"



TK-6-5202
Gardner
16 cm, 6 1/4"




TK-6-5203
Liston-mini
15 cm, 6"

Page 128–129


- Bone Elevators, Bone Files and Chisels





TK-6-5204
Williger
16 cm, 6 1/4"

TK-6-5205
Henahan
17 cm, 6 3/4"




TK-6-5206
Freer
19 cm, 7 1/2"

TK-6-5207
Bennet
17 cm, 6 3/4"

TK-6-5208
Seldin
19 cm, 7 1/2"

TK-6-5209
Seldin
19,5 cm, 7 3/4"

Bone Files and Chisels



Page 130–151

- Pliers For Orthodontics
- Mouth Gags
- Tongue Depressors
- Lip and Cheek Retractors
- Implantology Instrumetns



Pinze per formare loops
Loops forming pliers



TK-7-5101
Adams

Adams pliers for bending
wires up to 0,7mm. (.28").



TK-7-5102
Adams

Adams pliers for bending
wires up to 0,7mm. (.28").



TK-7-5103
Adams

Tweed loop forming pliers.
For forming loop in wires up
to 0,56 mm. (.22").



TK-7-5104
Adams

Aderer three jaws bending
pliers. For wires up
to 0,7 mm. (.28").



TK-7-5105
Adams

Aderer three jaws bending
pliers. For wires up to 0,7 mm.
(.28").

Ligature Cutter Pliers



TK-7-5106

Micro-miniature light wire cutter. It cuts soft ligature wires. For use in caso of limited access. Caution: maximum cutting capacity: soft ligature wire up to 0,30 mm. (.012").



TK-7-5107

Wire cutter. For ligature wires and pins. It is strong enough to bend archwire ends after distal cutting. Maximum cutting capacity up to 0,56 x 0,7 mm. (.022" x .028").



TK-7-5108

Wire cutter. For ligature wires and pins. It is strong enough to bend archwire ends after distal cutting. Maximum cutting capacity: up to 0,9 mm. (.036").



TK-7-5109

Wire cutter. For ligature wires and pins. It is strong enough to bend archwire ends after distal cutting. Maximum cutting capacity: up to 1,5 mm. (.044").

Pinze per formare laboratorio
Laboratory pliers



TK-7-5110

Flat nose pliers with smooth beaks for bending wires up to 0,7 mm. (.028").



TK-7-5111

Bending wires up to 0,7 mm. (.028").



TK-7-5112

Bending wires up to 0,7 mm. (.028").



TK-7-5113

Bending wires up to 0,7 mm. (.028").




TK-7-5114
Nance

Nance loop closing pliers, for bending of stop loops and special for Adams' clasps. Wires up to 0,7 mm. (.028").



TK-7-5115
Nance

Universal pliers for bending wires up to 0,9 mm. (.035") or cutting up to 0,7 mm. (.028").



TK-7-5116

Stainless steel pliers tray standard dimensions



TK-7-5117
Universal

Universal pin roughening, cutting and bending pliers. For wires up to 0,9 mm. (.035"). Will cut wires up to 0,7 mm. (.028").

Utility and specialty pliers



TK-7-5118
Reynolds

Reynolds contouring pliers.



TK-7-5119
Reynolds

Lingual arch pliers used to form double and triple back bends in 0,7 mm. (.028") or 0,9 mm. (.035") wire for insertion in lingual sheaths.



TK-7-5120
Weingart

Weingart pliers with serrated points.



TK-7-5121
De La Rosa

De La Rosa arch forming pliers. Designed to form and contour arch wires on round or rectangular wire. The beak has grooves of 0,4 mm (.016"), 0,46 mm (.018") and 0,56 mm (.022").



TK-7-5122
De La Rosa

Arch forming pliers with a non grooved contouring surface the forming of precise arches.

Loop formare pliers



TK-7-5123
Waldsachs

Waldsachs universal pliers for bending wires up to 1,0 mm. (.039") or cutting up to 0,7 mm. (.028").



TK-7-5124

Round and Concave pliers for loop forming and bending of thin round and square wires up to 0,5 mm. (.020").



TK-7-5125
Tweed

Tweed loop Forming Pliers. For round or rectangular wires up to 0,56 x 0,7 mm. (.022" x .028"). Replaceable cylindrical beak with three sizes: 1,15 mm. (.045), 1,50 mm. (.059"), 1,90 mm. (.075"). Concave beak with slight parallel serrations, to keep the wire in proper position.



TK-7-5126
Tweed

Tweed Loop Forming Pliers. For making perfect Omega Loops with round or rectangular wires up to 0,56 x 0,7 mm. (.022" x .028"). The concave beak is exactly the same as pliers 5809-1



TK-7-5127
Tweed

Single diameter loop forming pliers. Wire up to 0,56 x 0,7 mm. (.022" x 0.28"). Replaceable concave beak has only one diameter section of 1,50 mm. (.059").

Distal end cutters



TK-7-5128

Distal end cutter -safety hold- a precision cutter designed for distal arch cutting in the patient's mouth. It cuts hard wires end holds the loose distal end preventing any damaging of oral tissue. Maximum cutting capacity rectangular wires up to 0,56 x 0,7 mm. (.022" x .020").




TK-7-5129

Micro- distal end cutter - safety hold-. Similar to 3227



TK-7-5130

Direct bonding bracket remover (wide). The removal grip firmly for easy incisal-lingival (vertical) use or, if needed, for mesial-distal (horizontal) use with the same ease.



TK-7-5131


Direct bond attachment remover (angled). For removing all bonded attachments. May be used mesial-distal-horizontal) or incisal-lingival (vertical) in anterior areas. Angled perfectly for easy access to difficult posterior areas.

Bird beack pliers



TK-7-5132

Wire bending pliers. For wires up to 0,9 mm. (.036").



TK-7-5133

Bird Beak Angle pliers for bending round, square or rectangular wire up to 0,7 mm. (.028").



TK-7-5134

Double rounded jaw wire bending pliers. For wires up to 0,7 mm. (.028").



TK-7-5135

Bird beak pliers for bending round, square or rectangular wire up to 0,7 mm. (.028"). With cutter.



TK-7-5136


Double rounded jaw wire bending pliers. For wires up to 0,7 mm. (.028").



TK-7-5137
USA Model
125 mm




TK-7-5138
Waldron
125 mm



TK-7-5139
How
140 mm



TK-7-5140
How
145 mm




TK-7-5141
Adams
125 mm



TK-7-5142
Angle
125 mm



TK-7-5143
young
130 mm



TK-7-5144
Nance
130 mm




TK-7-5145
Weingart
140 mm




TK-7-5146
Johnson
130 mm




TK-7-5147
Mc Kellops
155 mm



TK-7-5148
Universal
150 mm



TK-7-5149
Waldsachs
155 mm



TK-7-5150
Aderer
125 mm



TK-7-5151
Nance
130 mm



TK-7-5152
Langenbeck
140 mm



TK-7-5153
Langenbeck
130 mm



TK-7-5154
Langenbeck
130 mm



TK-7-5155
Peeso
125 mm



TK-7-5156
Goslee
125 mm




TK-7-5157
Fischer
140 mm




TK-7-5158
Andresen
130 mm



TK-7-5159
Schwarz
130 mm



TK-7-5160
MSY
130 mm



TK-7-5161
Langenbeck
125 mm



TK-7-5162



TK-7-5163



TK-7-5164



TK-7-5165



TK-7-5166




TK-7-5167




TK-7-5168




TK-7-5169




TK-7-5170




TK-7-5171




TK-7-5172




TK-7-5173




TK-7-5174




TK-7-5175



TK-7-5176




TK-7-5177




TK-7-5178




TK-7-5179




TK-7-5180




TK-7-5181




TK-7-5182




TK-7-5183



TK-7-5184



TK-7-5185



TK-7-5186



Molt
TK-7-5187 10 cm
TK-7-5188 11 cm
TK-7-5189 14 cm



TK-7-5190 16 cm
TK-7-5191 19 cm
 Roser-König



TK-7-5192
 Denhardt
 13 cm



TK-7-5193
 Fergusson - Ackland
 14.5 cm



TK-7-5194
 Heister
 Maunier Wood
 9.5 cm




TK-7-5195
 Obwegeser 11 cm
 Spreading 47 mm



TK-7-5196
 Heister
 13 cm



TK-7-5197 Small
TK-7-5198 Medium
TK-7-5199 Large
 Cheek Retractor



TK-7-5200
Hartmann
15 cm



TK-7-5201
Doyen Tongue Spatula
13.5 cm




TK-7-5202
14.5 cm




TK-7-5203
18 cm




TK-7-5204
Bruenings
19 cm




TK-7-5205
Tobold
15 cm




TK-7-5206
Tobold
16 cm



TK-7-5207
14 cm



TK-7-5208
15 cm



Lip and Cheek Retractors



TK-7-5209
Obwegeser Chin
Retractor
16.5 cm



Obwegeser 16 cm
TK-7-5210 8 mm
TK-7-5211 10 mm



TK-7-5212 43 x 11 mm
TK-7-5213 57 x 11 mm
TK-7-5214 83 x 16 mm

Obwegeser soft
Tissue Retractor
21.5 cm



TK-7-5215
Obwegeser 35 x 11 mm
Retractor
22 cm



TK-7-5216
Obwegeser
Ramus retractor for
separation of ramus & tissue
22 cm



TK-7-5217
Simplex
10.5 cm



30x19 mm

33x25 mm

TK-7-5218
Kilner
15 cm



TK-7-5219
Minnesota
14 cm



TK-7-5220
Cawood-Minnesota
15 cm



TK-7-5221
Austin



TK-7-5222
Hilger
Fits any simple stem mouth mirror handle



TK-7-5223
Hilger
with handle




TK-7-5224
Max
20 cm

Finger protectors



TK-7-5225 Small
TK-7-5226 Medium
TK-7-5227 Large

Langenbeck



TK-7-5228
Oringer
11.5 cm



TK-7-5229
Obwegeser Cheek Retractor
12cm



TK-7-5230
Columbia
14.5 cm



Sternberg


TK-7-5231 12 cm
TK-7-5232 14 cm
TK-7-5233 15 cm



Vestibulum
TK-7-5234




TK-7-5235
Max



TK-7-5236
Max

Implants Retractors

Retractors



TK-7-5237

TK-7-5238

Senn-Miller
16 cm



TK-7-5239

TK-7-5240

TK-7-5241

TK-7-5242

Max
16 cm




TK-7-5243

TK-7-5244


TK-7-5245

TK-7-5246

Max
16 cm



Kirkland
Gum & Tissue
Retractors



TK-7-5247
Fig. 23 K



TK-7-5248
Fig. 24 K



TK-7-5249
Fig. 27 K



TK-7-5250
Fig. 28 K



TK-7-5251
Collin Tongue holding forceps
17 cm



TK-7-5254
Rubber Pad pair

Young Tongue holding forceps
TK-7-5252 15 cm
TK-7-5253 17 cm



TK-7-5255
Farabeuf 12 cm
Pair



TK-7-5256
Parker-Langenbeck 21 cm
Pair

Osteotomes

# 1	Diam.	# 2	Diam.	# 3	Diam.	# 4	Diam.	# 5	Diam.	# 6	Diam.
0	Ø 1.60	0	Ø 1.90	0	Ø 2.80	0	Ø 3.10	0	Ø 3.90	0	Ø 5.50
7.5	Ø 1.80	7.5	Ø 2.20	7.5	Ø 3.20	7.5	Ø 3.70	10	Ø 5.00	10	Ø 6.20
10	Ø 2.00	10	Ø 2.35	10.5	Ø 3.35	10.5	Ø 3.90	13	Ø 5.20	13	Ø 6.30
13.5	Ø 2.40	13.5	Ø 3.30	13.5	Ø 3.60	15.5	Ø 4.30	20.5	Ø 6.30		
15.5	Ø 3.60	15.5	Ø 3.40	21	Ø 4.40						
		23.5	Ø 4.50								



15.5 cm

1
Ø 1.6 mm



TK-7-5257
15.5 cm

2
Ø 1.9 mm



TK-7-5258
15.5 cm

3
Ø 2.8 mm



TK-7-5259
15.5 cm

4
Ø 3.10 mm



TK-7-5260
15.5 cm

5
Ø 3.90 mm



TK-7-5261
15.5 cm

6
Ø 5.50 mm



TK-7-5262
15.5 cm

# 1	Diam.	# 2	Diam.	# 3	Diam.	# 4	Diam.	# 5	Diam.
0	Ø 1.50	0	Ø 2.10	0	Ø 2.60	0	Ø 3.00	0	Ø 3.50
5	Ø 1.90	5	Ø 2.60	5	Ø 3.00	5	Ø 3.50	5	Ø 3.90
8	Ø 2.20	8	Ø 2.80	8	Ø 3.30	8	Ø 3.80	8	Ø 4.20
10	Ø 2.40	10	Ø 3.00	10	Ø 3.50	10	Ø 4.00	10	Ø 4.40
12	Ø 2.60	12	Ø 3.20	12	Ø 3.70	12	Ø 4.20	12	Ø 4.60
14	Ø 2.80	14	Ø 3.40	14	Ø 3.90	14	Ø 4.40	14	Ø 4.80

1
Ø 1.5 mm



TK-7-5263
15.5 cm

2
Ø 2.1 mm



TK-7-5264
15.5 cm

3
Ø 2.6 mm



TK-7-5265
15.5 cm

4
Ø 3.0 mm



TK-7-5266
15.5 cm

5
Ø 3.5 mm



TK-7-5267
15.5 cm

Ø 2.8 mm



TK-7-5268
15.5 cm

Ø 3.5 mm



TK-7-5269
15.5 cm


Ø 4.2 mm




TK-7-5270
15.5 cm

Spelling Straight
Bone Pusher

Osteotomes



Convex Tips



Osteotomes with Stop

Concave

Ø2mm Ø3mm Ø3.3mm Ø4mm

6mm 11mm 15mm

TK-7-5292 TK-7-5293 TK-7-5294 TK-7-5295

Convex

Ø2mm Ø3mm Ø3.3mm Ø4mm

6mm 11mm 15mm

TK-7-5296 TK-7-5297 TK-7-5298 TK-7-5299

Pointed

Ø2mm Ø3mm Ø3.3mm Ø4mm

6mm 11mm 15mm

TK-7-5300 TK-7-5301 TK-7-5302 TK-7-5303

16.5 cm Str

Concave

Ø2.5mm Ø3mm Ø3.8mm Ø4.5mm Ø5mm

8mm 11mm 13mm 15mm 17mm 20mm

TK-7-5316 TK-7-5317 TK-7-5318 TK-7-5319 TK-7-5320

Convex

Ø2.5mm Ø3mm Ø3.8mm Ø4.5mm Ø5mm

8mm 11mm 13mm 15mm 17mm 20mm

TK-7-5321 TK-7-5322 TK-7-5323 TK-7-5324 TK-7-5325

Pointed

Ø2.5mm Ø3mm Ø3.8mm Ø4.5mm Ø5mm

8mm 11mm 13mm 15mm 17mm 20mm

TK-7-5326 TK-7-5327 TK-7-5328 TK-7-5329

16.5 cm Str

Osteotomes with Stop

Concave

Ø2mm Ø3mm Ø3.3mm Ø4mm

6mm 11mm 15mm

TK-7-5304 TK-7-5305 TK-7-5306 TK-7-5307

Convex

Ø2mm Ø3mm Ø3.3mm Ø4mm

6mm 11mm 15mm

TK-7-5308 TK-7-5309 TK-7-5310 TK-7-5311

Pointed

Ø2mm Ø3mm Ø3.3mm Ø4mm

6mm 11mm 15mm

TK-7-5312 TK-7-5313 TK-7-5314 TK-7-5315

16.5 cm Cvd

Concave

Ø2.5mm Ø3mm Ø3.8mm Ø4.5mm Ø5mm

8mm 11mm 13mm 15mm 17mm 20mm

TK-7-5330 TK-7-5331 TK-7-5332 TK-7-5333 TK-7-5334

Convex

Ø2.5mm Ø3mm Ø3.8mm Ø4.5mm Ø5mm

8mm 11mm 13mm 15mm 17mm 20mm

TK-7-5335 TK-7-5336 TK-7-5337 TK-7-5338 TK-7-5339

Pointed

Ø2.5mm Ø3mm Ø3.8mm Ø4.5mm Ø5mm

8mm 11mm 13mm 15mm 17mm 20mm

TK-7-5340 TK-7-5341 TK-7-5342 TK-7-5343 TK-7-5344

16.5 cm Cvd

PDL Elevators



TK-7-5345 Distal, 2.5 mm Micro Edge



TK-7-5346 Mesial, 2.5 mm Micro Edge



TK-7-5347 Distal, 4 mm Micro Edge



TK-7-5348 Mesial 4 mm Micro Edge



TK-7-5349 2.5 mm Micro Edge, Straight



TK-7-5350 2.5 mm Micro Edge, Curved



TK-7-5351 3.4 mm

Straight **TK-7-5352** Set

Of 7 Pieces



TK-7-5353 Distal, 2.5 mm Micro Edge



TK-7-5354 Mesial, 2.5 mm Micro Edge




TK-7-5355 Distal, 4 mm Micro Edge



TK-7-5356 Mesial 4 mm Micro Edge



TK-7-5357 2.5 mm Micro Edge, Straight



TK-7-5358 2.5 mm Micro Edge, Curved



TK-7-5359 3.4 mm Straight

TK-7-5360 Set Of 7 Pieces



Instruments Tray

Trephine Drills

							
TK-7-5361 Interior 2 mm Exterior 3 mm	TK-7-5362 Interior 3 mm Exterior 4 mm	TK-7-5363 Interior 4 mm Exterior 5 mm	TK-7-5364 Interior 5 mm Exterior 6 mm	TK-7-5365 Interior 6 mm Exterior 7 mm	TK-7-5366 Interior 7 mm Exterior 8 mm	TK-7-5367 Interior 8 mm Exterior 9 mm	TK-7-5368 Interior 9 mm Exterior 10 mm




TK-7-5369
Trephine Drills Stand
Complete with drills


Tissue Punches



TK-7-5370
Trephine Drills Stand
Complete with drills

		
TK-7-5371 2 mm	TK-7-5375 6 mm	TK-7-5378 3 mm
		
TK-7-5372 3 mm	TK-7-5376 7 mm	TK-7-5379 4 mm
		
TK-7-5373 4 mm	TK-7-5377 8 mm	TK-7-5380 5 mm
		
TK-7-5374 5 mm		TK-7-5381 6 mm

Keyes



TK-7-5382
Bone Collector Chisel Straight



TK-7-5383
Blade for Straight bone collector chisel



TK-7-5384
Bone Collector Chisel Curved



TK-7-5385
Blade for Curved bone collector chisel



TK-7-5386
Ratchet Screw driver handle
Hexagonal head



TK-7-5387
Bone Morselizer with Titanium Grinding



TK-7-5388
Bone Mixing Basin
Dia 35 mm



TK-7-5389
Manual Easy Bone Mill



TK-7-5390
Bone Craker Titanium

Page 152–155

- Impression Trays



All Impression Trays Upper & Lower Type
are available in S, M, L & XL Sizes



TK-7-5391
"HYDRO-LOCK"
stainless steel
water-cooled
impression trays
with retention rim.



TK-7-5392
"HYDRO-LOCK"
stainless steel
water-cooled impression
trays with retention rim.



TK-7-5393
"NEW SUPER"
stainless steel
impression trays
with retention rim.



TK-7-5394
"NEW SUPER" stainless
steel impression trays
with retention rim.



TK-7-5395
"ANATOMIC" solid
regular stainless
steel impression
trays without
retention rim.



TK-7-5396
"ANATOMIC" solid regular
stainless steel impression
trays without retention rim.



TK-7-5397
Perforated-total denture
stainless steel
impression trays with rim.



TK-7-5398
Perforated-total denture
stainless steel
impression trays with rim.

Stainless Steel Impression Trays

All Impression Trays Upper & Lower Type
are available in S, M, L & XL Sizes



TK-7-5399
"PERMA-LOCK" solid regular stainless
steel impression trays with retention rim.



TK-7-5400
"PERMA-LOCK" solid regular stainless
steel impression trays with retention rim.



TK-7-5401



TK-7-5402



TK-7-5403



TK-7-5404



TK-7-5405



TK-7-5406



TK-7-5407
Set of six

PACKINGS



TK-7-5408
Set of Eight



TK-7-5409
Solid stainless steel
anterior depressed
impression trays with
retention rim.



TK-7-5410



TK-7-5411



TK-7-5412



TK-7-5413



TK-7-5414



TK-7-5415



TK-7-5416
Solid stainless steel
anterior depressed
impression trays with
retention rim.

All Impression Trays Upper & Lower Type
are available in S, M, L & XL Sizes



TK-7-5417
Fig. 0 per bambini - for children



TK-7-5418
For narrow toothed upper jaws



TK-7-5419
Fig. 0 per bambini - for children



TK-7-5420
For partially toothed upper jaws



TK-7-5421
For partially toothed lower jaws



TK-7-5422
For edentulous upper jaws



TK-7-5423
For edentulous lower jaws

Stainless Steel Impression Trays

All Impression Trays Upper & Lower Type
are available in S, M, L & XL Sizes



TK-7-5424
Solid total denture stainless steel impression
tray with rim.



TK-7-5425
Solid total denture stainless steel impression
tray with rim.



TK-7-5426
"PERMA-LOCK" perforated regular stainless
steel impression trays with retention rim.



TK-7-5427
"PERMA-LOCK" perforated regular stainless
steel impression trays with retention rim.



TK-7-5428
Perforated stainless steel anterior depressed
impression trays with retention rim.



TK-7-5429



TK-7-5430



TK-7-5431



TK-7-5432



TK-7-5433



TK-7-5434



TK-7-5435
Perforated stainless steel anterior depressed
impression trays with retention rim.



TK-7-5436
Perforated stainless
steel impression
trays without
retention rim for
complete denture.



TK-7-5437
Perforated stainless steel
impression
trays without
retention
rim for complete denture.

Page 156–157

- Arriculators





TK-8-5101




TK-8-5102



TK-8-5103



TK-8-5104



TK-8-5105



TK-8-5106



TK-8-5107



TK-8-5108



TK-8-5109




TK-8-5110

Page 158–160

- Dental Surgical Hollowares





TK-8-5111
ENDODONTICBOX



TK-8-5112
STERILIZINGBOX
Available in dif. sizes
17x17x3 cm
20x10x4 cm
25x12x6 cm
30x14x6 cm



TK-8-5113 ENDO
TRAY BOX



TK-8-5114
COTTONDISPENSER



TK-8-5115
ALCOHOLLAMP

Sterilizing Drums



TK-8-5116



TK-8-5117



TK-8-5118

Sterilizing Boxes



TK-8-5119



TK-8-5120



TK-8-5121

Sterilizing Boxes



TK-8-5122



TK-8-5123



TK-8-5124



TK-8-5125



TK-8-5126



TK-8-5127

FORCEP AND THERMOMETER JARS

COTTON DISPENSER



TK-8-5128
IRRIGATOR(Douchcan)



TK-8-5129
MEASURINGJUG



TK-8-5130
KIDNEY TRAY



TK-8-5131
FORCEP JAR



TK-8-5132
GALLIPOT WITH & WITHOUT LID



TK-8-5133
SPITTEN MUG



TK-8-5134
SPIRIT LAMP



TK-8-5135
Perforated Tray



TK-8-5136
Perforated Base



TK-8-5137
Solid Base



TK-8-5138
Instruments Tray



TK-8-5139
Instruments Tray



TK-8-5140
**Instruments Tray
(Perforated)**



TK-8-5141
**Instruments Tray Shallow
(Mayo)**



TK-8-5142
Scalers Tray



TK-8-5143
Cathetar Tray

Alphabetical Index

A					
Adams	125, 133	Cleveland	122	Goslee	134
Aderer	133	College	105, 107	Gouges	68
Adson	96	Columbia	64, 66, 74	Gracey Curettes	69
Ainsworth	110	Columbia	79	Gracey Mesial Distal	69
Albee	77	Cone-Socket	90	Gracey	56-69-72
Albrecht	109	Cork Screw Tip	104	Gritman	76-78-80
Alcohol Lamp	142	Cotton Dispenser	142, 143		
Alginate	93-94	Coupland	37	H	
Allen	109	Crane-kaplan	74	Halstead-Mosquito	118
Amalgam Carriers	117	Crile-Wood	96	Hegar-Baumgartner	95-96
Amalgam Guns	117	Crown Removers	115-116	Heidbrink	45-48-50
Amalgam Well	117	Cryer	35, 37, 48	Heidemann	55
Andresen	134	Curette	64, 66	Hemingway	75
Angle	133			Henahan	123
Antkoviak	78	D		How	133
Apical	36, 44	De La Rosa	129	Hylin	40, 76
Arruga	97	Dental Probes	57		
Articulators	140	Dental Syringes	103-104	I	
Ash	53, 88	Derf	95	Impression Trays	136-139
Aspirating Syringes	102	Didner	113	Intra-Ligamental Syringes	104
				Irrigator	143
B		E		Ivory	110, 113
Band Adaptor	102	Eliselberg-Mathieu	97	Iwanson	77
Barry	47, 50	Elliott	53		
Baumgarthner	96	Endo Finger Ruler	102	J	
Beale	76	Endo Tray Box	142	Jaquette	74
Beebee	116	Endodontic Box	142	Johnson	133
Bein	33, 40	Ermert	78		
Bennet	53, 123	Excavators	87-88	K	
Beyer	121	Explorers	81	Kelly	118
Bicuspids	114	Extracting Forceps	01-31	Kerrison	120
Birdbeak Pliers	132			Kidney Tray	143
Blumenthal	120, 121	F		Kilian	124
Bohler	120	Fahnenstock	79	Kilian-Claus	124
Bohm	118	Filling Instruments	52-56	Klein	07
Boley	77	Fischer	134	Kleinert-Kutz	121
Bone Chisels	68	Flagg	107	Kopp	39, 49, 108
Boynton	97	Flohr	38		
Breech Loading Syringe	102	Forcepps Jar	143, 144	L	
Brewer	110	Fox	58	Landgraf	107
Buck	61-63, 65	Fraizer	100	Langenbeck	134
Buckley	67, 124	Freer	123	Le Cron	76
Burnishers	58	Friedman	119-120	Lecluse	43
		Furrer	119	Lessmann	79
C		G		Lichtenberg	97
Caroll	77	Gallipot	144	Listion-mini	122
Castroviejo	77, 98, 100	Gardner	67, 122	London Hospital	41
Cathetar Tray	144	GF1	108	London-College	107, 109
Cement Spatulas	91, 92	GF2	108	Lucas	75
Chompret	44	Gillies	96	Luer	120
Ck1	108	Goldman Fox	60, 63, 65, 73, 120	Luks	84
CK2	108				

M

Machtou	83
Mallets	68
Mathieu	99, 100
Matrix Bands	113, 114
Matrix Retainers	113
Mayo-Hegar	95
Mc Call	73
Mc Kellops	133
Mead	43, 58, 119
Measuring Jug	143
Measuring Plate	84
Medicon	78
Meriam	106, 109
Merrit A	60
Mershon	110, 116
Micro-Barraquer	99-100
Micro-Mosquito	118
Miller	118, 124
Mini-Ryder	95
Molar	114
Morell	116
Mouth Mirror Handles	89, 90
Mouth Mirrors	89, 90
MSY	76, 134

N

Nabers	83
Nance	128, 133
Niro	122

O

Obwegeser	124
Ochsenbein	61, 62, 67
Olsen-Heger	96
Oms	83
Ott	41

P

P.K. Thomas	70
Partsch	68
Pocket Probes	134
Plaster Knives	58-60
Plaster Spatulas	80
Position Gauge	93-94
Pott	102
	42

R

Rein	81
Relax	05
Reynolds	129
Rhodes	61, 62

Rhodium

Rhodium	89
Roach	76
Root Elevators	32--51
Root Splinter	31
Root-Canal Excavators	86
Root-Canal Explorers	86
Root-Canal Plungers	84-85
Rubber Dam Clamps	111-112
Rubber Dam Frame	110
Ruskin-Mini	121

S

Scalers	58
Schell	119
Schilder	85
Schluger	61, 62
Schwarz	134
Seldin	34, 48, 123
Spirit Lamp	144
Spitten Mug	144
Spreader	82-83
Sterilizing Boxes	142, 143
Sterilizing Drums	142
Sugarman	61, 162

T

Taylor	74
Toennis	97
Tofflemire	114
Tofflemire-Junior	113
Tofflemire-Senior	113
Tofflemire-Universal	113
Towner Jaquette	74
Towner	74
Trays	144
Treymann	116
Tweed	130

U

Univarsal “HM”	116
Univarsal	128, 133
University of Michigon	60
USA Model	133

V

Vehe	76
Vernier	77
Vertical Condenser	83

W

Wakefield	
Waldron	67
Waldsachs	133
Warwick-James	130, 133
Water Syringe	41
Weingart	101
White	129, 133
Williams	52, 60, 83
Williger	123
Winter	46, 50
Woodson	54, 55

Y

Young	133
Younger-Good	73

Z

Zahle	76
Zaufal-Jansen	121
Zerfing-Chisel	64, 66
Zuricher	77

WE HAVE ALL RANGE OF INSTRUMENTS



CONTACT INFO



+92 344 6346401



info@tkmed-instruments.com



www.tkmed-instruments.com



Head Marala Road Machi Khokhar, Sialkot, 51310

Get your high quality Instruments